

Bovine AQP3/Aquaporin-3 ELISA Kit

Catalog No: #EK19318

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine AQP3/Aquaporin-3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Aquaporin-3
Target Name	Aquaporin-3
Other Names	AQP-3, Aquaglyceroporin-3
Accession No.	Q08DE6
Uniprot	Q08DE6
GeneID	780866;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/ml Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: AQP3

Note: This product is for in vitro research use only