

Pig AQP3/Aquaporin-3 ELISA Kit

Catalog No: #EK19319

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pig AQP3/Aquaporin-3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Aquaporin-3
Target Name	Aquaporin-3
Other Names	Aquaglyceroporin-3
Accession No.	A9Y006
Uniprot	A9Y006
GeneID	100126235;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/ml Sensitivity:0.088ng/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:AQP3

Note: This product is for in vitro research use only