

Rat Aqp3/Aquaporin-3 ELISA Kit

Catalog No: #EK19320

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Aqp3/Aquaporin-3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Aquaporin-3
Target Name	Aquaporin-3
Other Names	AQP-3,31.4 kDa water channel protein,Aquaglyceroporin-3
Accession No.	P47862
Uniprot	P47862
GeneID	65133;

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.19ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Aqp3

Note: This product is for in vitro research use only