

Bovine AQP4/Aquaporin-4 ELISA Kit

Catalog No: #EK19321

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine AQP4/Aquaporin-4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Aquaporin-4
Target Name	Aquaporin-4
Other Names	AQP-4,Mercurial-insensitive water channel,MIWC,WCH4
Accession No.	O77750
Uniprot	O77750
GeneID	281008;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.161ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:AQP4

Note: This product is for in vitro research use only