

## Pig AQP5/Aquaporin-5 ELISA Kit

Catalog No: #EK19322

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pig AQP5/Aquaporin-5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Aquaporin-5
Target Name	Aquaporin-5
Other Names	AQP-5
Accession No.	A8W649
Uniprot	A8W649
GeneID	100126278;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.052ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:AQP5

Note: This product is for in vitro research use only