

Rat Aqp5/Aquaporin-5 ELISA Kit

Catalog No: #EK19323

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Aqp5/Aquaporin-5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Aquaporin-5
Target Name	Aquaporin-5
Other Names	AQP-5
Accession No.	P47864
Uniprot	P47864
GeneID	25241;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.164ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Aqp5

Note: This product is for in vitro research use only