

## Bovine CCND1/G1/S-specific cyclin-D1 ELISA Kit

Catalog No: #EK19324

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine CCND1/G1/S-specific cyclin-D1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine G1/S-specific cyclin-D1
Target Name	G1/S-specific cyclin-D1
Accession No.	Q2KI22
Uniprot	Q2KI22
GeneID	524530;

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.165ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CCND1

Note: This product is for in vitro research use only