

Mouse Ccnd2/G1/S-specific cyclin-D2 ELISA Kit

Catalog No: #EK19325

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ccnd2/G1/S-specific cyclin-D2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse G1/S-specific cyclin-D2
Target Name	G1/S-specific cyclin-D2
Other Names	Cyl-2
Accession No.	P30280
Uniprot	P30280
GeneID	12444;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.081ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ccnd2

Note: This product is for in vitro research use only