

# Mouse Akr1c6/Estradiol 17 beta-dehydrogenase 5 ELISA Kit

Catalog No: #EK19328

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Akr1c6/Estradiol 17 beta-dehydrogenase 5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Estradiol 17 beta-dehydrogenase 5
Target Name	Estradiol 17 beta-dehydrogenase 5
Other Names	17-beta-HSD 5,Hsd17b5
Accession No.	P70694
Uniprot	P70694
GeneID	83702;

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39.4pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Akr1c6

Note: This product is for in vitro research use only