

Rat Cpb2/Carboxypeptidase B2 ELISA Kit

Catalog No: #EK19329

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Rat Cpb2/Carboxypeptidase B2 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Rt |
| Specificity | Natural and recombinant Rat Carboxypeptidase B2 |
| Target Name | Carboxypeptidase B2 |
| Other Names | Carboxypeptidase R,CPR,Carboxypeptidase U,CPU,Thrombin-activable fibrinolysis inhibitor,TAFI |
| Accession No. | Q9EQV9 |
| Uniprot | Q9EQV9 |
| GeneID | 113936; |

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cpb2

Note: This product is for in vitro research use only