

Human ALDOB/Fructose-bisphosphate aldolase B ELISA Kit



Catalog No: #EK19330

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ALDOB/Fructose-bisphosphate aldolase B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Fructose-bisphosphate aldolase B
Target Name	Fructose-bisphosphate aldolase B
Other Names	Liver-type aldolase,ALDB
Accession No.	P05062
Uniprot	P05062
GeneID	229;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.5pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ALDOB

Note: This product is for in vitro research use only