

Human CASP9/Caspase-9 ELISA Kit

Catalog No: #EK19333

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CASP9/Caspase-9 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Caspase-9
Target Name	Caspase-9
Other Names	CASP-9, Apoptotic protease Mch-6, Apoptotic protease-activating factor 3, APAF-3, ICE-like apoptotic protease 6, ICE-LAP6, MCH6
Accession No.	P55211
Uniprot	P55211
GeneID	842;

Application Details

Assay Type: Sandwich Detect Range: 3.12-200ng/ml Sensitivity: 1.64ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CASP9

Note: This product is for in vitro research use only