

# Mouse Bid/BH3-interacting domain death agonist ELISA Kit



Catalog No: #EK19334

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Bid/BH3-interacting domain death agonist ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse BH3-interacting domain death agonist
Target Name	BH3-interacting domain death agonist
Other Names	p22 BID,BID
Accession No.	P70444
Uniprot	P70444
GeneID	12122;

## Application Details

Assay Type: Sandwich Detect Range:1.56-100ng/mL Sensitivity:0.78ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Bid

Note: This product is for in vitro research use only