

Mouse Cfh/Complement factor H ELISA Kit

Catalog No: #EK19335

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cfh/Complement factor H ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Complement factor H
Target Name	Complement factor H
Other Names	Protein beta-1-H,Hf1
Accession No.	P06909
Uniprot	P06909
GeneID	12628;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.083ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cfh

Note: This product is for in vitro research use only