

Mouse Calm1/Calmodulin ELISA Kit

Catalog No: #EK19336

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Calm1/Calmodulin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Calmodulin
Target Name	Calmodulin
Other Names	CaM, Calm, Cam, Cam1
Accession No.	P62204
Uniprot	P62204

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.084ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Calm1

Note: This product is for in vitro research use only