

## Mouse Bmp6/Bone morphogenetic protein 6 ELISA Kit

Catalog No: #EK19337

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Bmp6/Bone morphogenetic protein 6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Bone morphogenetic protein 6
Target Name	Bone morphogenetic protein 6
Other Names	BMP-6, VG-1-related protein, VGR-1, Bmp-6, Vgr1
Accession No.	P20722
Uniprot	P20722
GeneID	12161;

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/ml Sensitivity: 0.33ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Bmp6

Note: This product is for in vitro research use only