

Rat Bmp6/Bone morphogenetic protein 6 ELISA Kit

Catalog No: #EK19338

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Bmp6/Bone morphogenetic protein 6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Bone morphogenetic protein 6
Target Name	Bone morphogenetic protein 6
Other Names	BMP-6,VG-1-related protein,VGR-1,Bmp-6,Vgr
Accession No.	Q04906
Uniprot	Q04906
GeneID	25644;

Application Details

Assay Type: Sandwich Detect Range:0.625-40ng/mL Sensitivity:0.325ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Bmp6

Note: This product is for in vitro research use only