

Human CD28/T-cell-specific surface glycoprotein CD28 ELISA Kit

Catalog No: #EK19339

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CD28/T-cell-specific surface glycoprotein CD28 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human T-cell-specific surface glycoprotein CD28
Target Name	T-cell-specific surface glycoprotein CD28
Other Names	TP44,CD28
Accession No.	P10747
Uniprot	P10747
GeneID	940;

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.445ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CD28

Note: This product is for in vitro research use only