

Mouse Cd28/T-cell-specific surface glycoprotein CD28 ELISA Kit



Catalog No: #EK19340

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cd28/T-cell-specific surface glycoprotein CD28 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse T-cell-specific surface glycoprotein CD28
Target Name	T-cell-specific surface glycoprotein CD28
Other Names	CD28
Accession No.	P31041
Uniprot	P31041
GeneID	12487;

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.446ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cd28

Note: This product is for in vitro research use only