

## Rat Aldh1a1/Retinal dehydrogenase 1 ELISA Kit

Catalog No: #EK19350

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Aldh1a1/Retinal dehydrogenase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Retinal dehydrogenase 1
Target Name	Retinal dehydrogenase 1
Other Names	RALDH 1,ALDH-E1,ALHDII,Aldehyde dehydrogenase family 1 member A1,Aldehyde dehydrogenase, cytosolic,Aldh
Accession No.	P51647
Uniprot	P51647
GeneID	24188;

## Application Details

Assay Type: Sandwich Detect Range:0.625-40ng/ml Sensitivity:0.328ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Aldh1a1

Note: This product is for in vitro research use only