

## Mouse Adh1/Alcohol dehydrogenase 1 ELISA Kit

Catalog No: #EK19351

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Mouse Adh1/Alcohol dehydrogenase 1 ELISA Kit          |
| Brief Description  | ELISA Kits  |
| Applications       | ELISA   |
| Species Reactivity | Ms  |
| Specificity        | Natural and recombinant Mouse Alcohol dehydrogenase 1 |
| Target Name        | Alcohol dehydrogenase 1                               |
| Other Names        | ADH-A2,Alcohol dehydrogenase A subunit,Adh-1          |
| Accession No.      | P00329  |
| Uniprot            | P00329  |
| GeneID             | 11522;  |

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39.1pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Adh1

Note: This product is for in vitro research use only