

Mouse Cd14/Monocyte differentiation antigen CD14 ELISA Kit



Catalog No: #EK19352

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cd14/Monocyte differentiation antigen CD14 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Monocyte differentiation antigen CD14
Target Name	Monocyte differentiation antigen CD14
Other Names	Myeloid cell-specific leucine-rich glycoprotein,CD14
Accession No.	P10810
Uniprot	P10810
GeneID	12475;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/ml Sensitivity:0.098ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name: Cd14

Note: This product is for in vitro research use only