

## Rat Calb2/Calretinin ELISA Kit

Catalog No: #EK19353

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Calb2/Calretinin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Calretinin
Target Name	Calretinin
Other Names	CR
Accession No.	P47728
Uniprot	P47728
GeneID	117059;

## Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 7.98pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Calb2

Note: This product is for in vitro research use only