

## Mouse Ccl24/C-C motif chemokine 24 ELISA Kit

Catalog No: #EK19355

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Ccl24/C-C motif chemokine 24 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse C-C motif chemokine 24
Target Name	C-C motif chemokine 24
Other Names	Eosinophil chemotactic protein 2, Eotaxin-2, Small-inducible cytokine A24, Scya24
Accession No.	Q9JKC0
Uniprot	Q9JKC0
GeneID	56221;

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/ml Sensitivity: 0.085ng/ml Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Ccl24

Note: This product is for in vitro research use only