

Mouse Apoe/Apolipoprotein E ELISA Kit

Catalog No: #EK19356

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Apoe/Apolipoprotein E ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Apolipoprotein E
Target Name	Apolipoprotein E
Other Names	Apo-E
Accession No.	P08226
Uniprot	P08226
GeneID	11816;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.064ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Apoe

Note: This product is for in vitro research use only