

Bovine CEBPA/CCAAT/enhancer-binding protein alpha ELISA Kit

Catalog No: #EK19357

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine CEBPA/CCAAT/enhancer-binding protein alpha ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine CCAAT/enhancer-binding protein alpha
Target Name	CCAAT/enhancer-binding protein alpha
Other Names	C/EBP alpha
Accession No.	O02754
Uniprot	O02754

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/ml Sensitivity: 0.084ng/ml Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CEBPA

Note: This product is for in vitro research use only