

# Mouse Cebpa/CCAAT/enhancer-binding protein alpha ELISA Kit



Catalog No: #EK19358

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cebpa/CCAAT/enhancer-binding protein alpha ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse CCAAT/enhancer-binding protein alpha
Target Name	CCAAT/enhancer-binding protein alpha
Other Names	C/EBP alpha,Cebp
Accession No.	P53566
Uniprot	P53566
GeneID	12606;

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39.5pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Cebpa

Note: This product is for in vitro research use only