

Bovine CEBPB/CCAAT/enhancer-binding protein beta ELISA Kit



Catalog No: #EK19359

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine CEBPB/CCAAT/enhancer-binding protein beta ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine CCAAT/enhancer-binding protein beta
Target Name	CCAAT/enhancer-binding protein beta
Other Names	C/EBP beta
Accession No.	O02755
Uniprot	O02755
GeneID	338319;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.064ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CEBPB

Note: This product is for in vitro research use only