

# Human CEBPB/CCAAT/enhancer-binding protein beta ELISA Kit



Catalog No: #EK19360

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CEBPB/CCAAT/enhancer-binding protein beta ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human CCAAT/enhancer-binding protein beta
Target Name	CCAAT/enhancer-binding protein beta
Other Names	C/EBP beta,Liver activator protein,LAP,Liver-enriched inhibitory protein,LIP,Nuclear factor NF-IL6,Transcription factor 5,TCF-5,TCF5,PP9092
Accession No.	P17676
Uniprot	P17676
GeneID	1051;

## Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.39ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CEBPB

Note: This product is for in vitro research use only