

Rat Cebpb/CCAAT/enhancer-binding protein beta ELISA Kit



Catalog No: #EK19361

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cebpb/CCAAT/enhancer-binding protein beta ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat CCAAT/enhancer-binding protein beta
Target Name	CCAAT/enhancer-binding protein beta
Other Names	C/EBP beta,C/EBP-related protein 2,Interleukin-6-dependent-binding protein,IL-6DBP,Liver-enriched inhibitory protein,LIP,Liver-enriched transcriptional activator,LAP,Silencer factor B,SF-B,Crp2,Nf-il6,Sfb
Accession No.	P21272
Uniprot	P21272
GeneID	24253;

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.23ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cebpb

Note: This product is for in vitro research use only