

Rat Bace1/Beta-secretase 1 ELISA Kit

Catalog No: #EK19362

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Bace1/Beta-secretase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Beta-secretase 1
Target Name	Beta-secretase 1
Other Names	Aspartyl protease 2,ASP2,Beta-site amyloid precursor protein cleaving enzyme 1,Beta-site APP cleaving enzyme 1,Memapsin-2,Membrane-associated aspartic protease 2,Bace
Accession No.	P56819
Uniprot	P56819
GeneID	29392;

Application Details

Assay Type: Sandwich Detect Range:1.56-100ng/mL Sensitivity:0.71ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Bace1

Note: This product is for in vitro research use only