

# Mouse Cebpe/CCAAT/enhancer-binding protein epsilon ELISA Kit

Catalog No: #EK19363

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cebpe/CCAAT/enhancer-binding protein epsilon ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse CCAAT/enhancer-binding protein epsilon
Target Name	CCAAT/enhancer-binding protein epsilon
Other Names	C/EBP epsilon,Gm294
Accession No.	Q6PZD9
Uniprot	Q6PZD9
GeneID	110794;

## Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.195ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Cebpe

Note: This product is for in vitro research use only