

Rat Cebpe/CCAAT/enhancer-binding protein epsilon ELISA Kit



Catalog No: #EK19364

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cebpe/CCAAT/enhancer-binding protein epsilon ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat CCAAT/enhancer-binding protein epsilon
Target Name	CCAAT/enhancer-binding protein epsilon
Other Names	C/EBP epsilon,C/EBP-related protein 1,Crp1
Accession No.	P56261
Uniprot	P56261
GeneID	25410;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.069ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cebpe

Note: This product is for in vitro research use only