

Rat Cdk5/Cyclin-dependent-like kinase 5 ELISA Kit

Catalog No: #EK19367

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cdk5/Cyclin-dependent-like kinase 5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cyclin-dependent-like kinase 5
Target Name	Cyclin-dependent-like kinase 5
Other Names	Cell division protein kinase 5, Serine/threonine-protein kinase PSSALRE, Tau protein kinase II catalytic subunit, TPKII catalytic subunit, Cdkn5
Accession No.	Q03114
Uniprot	Q03114
GeneID	140908;

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.393ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cdk5

Note: This product is for in vitro research use only