

# Rat C1qa/Complement C1q subcomponent subunit A ELISA Kit



Catalog No: #EK19368

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat C1qa/Complement C1q subcomponent subunit A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Complement C1q subcomponent subunit A
Target Name	Complement C1q subcomponent subunit A
Accession No.	P31720
Uniprot	P31720
GeneID	298566;

## Application Details

Assay Type: Sandwich Detect Range: 3.12-200ng/mL Sensitivity: 1.67ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: C1qa

Note: This product is for in vitro research use only