

Mouse Adam10/Disintegrin and metalloproteinase domain-containing protein 10 ELISA Kit



Catalog No: #EK19372

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Adam10/Disintegrin and metalloproteinase domain-containing protein 10 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Disintegrin and metalloproteinase domain-containing protein 10
Target Name	Disintegrin and metalloproteinase domain-containing protein 10
Other Names	ADAM 10,Kuzbanian protein homolog,Mammalian disintegrin-metalloprotease,CD156c,Kuz,Madm
Accession No.	O35598
Uniprot	O35598
GeneID	11487;

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.27ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Adam10

Note: This product is for in vitro research use only