

Rat Adam10/Disintegrin and metalloproteinase domain-containing protein 10 ELISA Kit

Catalog No: #EK19373

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Adam10/Disintegrin and metalloproteinase domain-containing protein 10 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Disintegrin and metalloproteinase domain-containing protein 10
Target Name	Disintegrin and metalloproteinase domain-containing protein 10
Other Names	ADAM 10,Kuzbanian protein homolog,Mammalian disintegrin-metalloprotease,CD156c,Madm
Accession No.	Q10743
Uniprot	Q10743
GeneID	29650;

Application Details

Assay Type: Sandwich Detect Range:1.56-100U/ml Sensitivity:1U/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Adam10

Note: This product is for in vitro research use only