

Chicken BCL2/Apoptosis regulator Bcl-2 ELISA Kit

Catalog No: #EK19374

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chicken BCL2/Apoptosis regulator Bcl-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken Apoptosis regulator Bcl-2
Target Name	Apoptosis regulator Bcl-2
Other Names	BCL-2
Accession No.	Q00709
Uniprot	Q00709
GeneID	396282;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.6pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: BCL2

Note: This product is for in vitro research use only