

Human CA2/Carbonic anhydrase 2 ELISA Kit

Catalog No: #EK19375

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CA2/Carbonic anhydrase 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Carbonic anhydrase 2
Target Name	Carbonic anhydrase 2
Other Names	Carbonate dehydratase II,Carbonic anhydrase C,CAC,Carbonic anhydrase II,CA-II
Accession No.	P00918
Uniprot	P00918
GeneID	760;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39.7pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CA2

Note: This product is for in vitro research use only