

Mouse Cdkn2a/Cyclin-dependent kinase inhibitor 2A, isoforms 1/2 ELISA Kit



Catalog No: #EK19378

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cdkn2a/Cyclin-dependent kinase inhibitor 2A, isoforms 1/2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cyclin-dependent kinase inhibitor 2A, isoforms 1/2
Target Name	Cyclin-dependent kinase inhibitor 2A, isoforms 1/2
Other Names	Cyclin-dependent kinase 4 inhibitor A,CDK4I,p16-INK4a,p16-INK4,P16ink4a
Accession No.	P51480
Uniprot	Q64364
GeneID	12578;

Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:8.06pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cdkn2a

Note: This product is for in vitro research use only