

Mouse Agt/Angiotensinogen ELISA Kit

Catalog No: #EK19379

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Agt/Angiotensinogen ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Angiotensinogen
Target Name	Angiotensinogen
Other Names	Serpin A8,Serpina8
Accession No.	P11859
Uniprot	P11859

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 8.08pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Agt

Note: This product is for in vitro research use only