

## Mouse Bmp7/Bone morphogenetic protein 7 ELISA Kit

Catalog No: #EK19381

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Bmp7/Bone morphogenetic protein 7 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Bone morphogenetic protein 7
Target Name	Bone morphogenetic protein 7
Other Names	BMP-7,Osteogenic protein 1,OP-1,Bmp-7,Op1
Accession No.	P23359
Uniprot	P23359
GeneID	12162;

## Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:4pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Bmp7

Note: This product is for in vitro research use only