

## Bovine CCK/Cholecystokinin ELISA Kit

Catalog No: #EK19382

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine CCK/Cholecystokinin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Cholecystokinin
Target Name	Cholecystokinin
Other Names	CCK
Accession No.	P41520
Uniprot	P41520
GeneID	617510;

## Application Details

Assay Type: Sandwich Detect Range: 6.25-400pg/mL Sensitivity: 3.235pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CCK

Note: This product is for in vitro research use only