

Dog CCK/Cholecystokinins ELISA Kit

Catalog No: #EK19383

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dog CCK/Cholecystokinins ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog Cholecystokinins
Target Name	Cholecystokinins
Other Names	CCK
Accession No.	Q9TS44
Uniprot	Q9TS44

Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 7.9pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CCK

Note: This product is for in vitro research use only