

Mouse Cck/Cholecystokinin ELISA Kit

Catalog No: #EK19384

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cck/Cholecystokinin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cholecystokinin
Target Name	Cholecystokinin
Other Names	CCK
Accession No.	P09240
Uniprot	P09240
GeneID	12424;

Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.7pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cck

Note: This product is for in vitro research use only