

## Pig CCK/Cholecystokinin ELISA Kit

Catalog No: #EK19385

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                             |
|--------------------|---------------------------------------------|
| Product Name       | Pig CCK/Cholecystokinin ELISA Kit           |
| Brief Description  | ELISA Kits                                  |
| Applications       | ELISA                                       |
| Species Reactivity | Pig                                         |
| Specificity        | Natural and recombinant Pig Cholecystokinin |
| Target Name        | Cholecystokinin                             |
| Other Names        | CCK                                         |
| Accession No.      | P01356                                      |
| Uniprot            | P01356                                      |
| GeneID             | 397468;                                     |

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ug/mL Sensitivity:0.156ug/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CCK

Note: This product is for in vitro research use only