

# Dog BEP/Beta-endorphin ELISA Kit

Catalog No: #EK19386



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Dog BEP/Beta-endorphin ELISA Kit           |
| Brief Description  | ELISA Kits                                 |
| Applications       | ELISA                                      |
| Species Reactivity | Dog  |
| Specificity        | Natural and recombinant Dog Beta-endorphin |
| Target Name        | Beta-endorphin                             |

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.085ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: BEP

Note: This product is for in vitro research use only