

Mouse Ctf1/Cardiotrophin-1 ELISA Kit

Catalog No: #EK19387

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ctf1/Cardiotrophin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cardiotrophin-1
Target Name	Cardiotrophin-1
Other Names	CT-1
Accession No.	Q60753
Uniprot	Q60753
GeneID	13019;

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.32ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ctf1

Note: This product is for in vitro research use only