

# Human ALOX15/Arachidonate 15-lipoxygenase ELISA Kit



Catalog No: #EK19389

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ALOX15/Arachidonate 15-lipoxygenase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Arachidonate 15-lipoxygenase
Target Name	Arachidonate 15-lipoxygenase
Other Names	15-LOX, 12/15-lipoxygenase, Arachidonate 12-lipoxygenase, leukocyte-type, 12-LOX, Arachidonate omega-6 lipoxygenase, LOG15
Accession No.	P16050
Uniprot	P16050
GeneID	246;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 35pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: ALOX15

Note: This product is for in vitro research use only