

## Bovine CST3/Cystatin-C ELISA Kit

Catalog No: #EK19390

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine CST3/Cystatin-C ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Cystatin-C
Target Name	Cystatin-C
Other Names	Colostrum thiol proteinase inhibitor, Cystatin-3
Accession No.	P01035
Uniprot	P01035
GeneID	281102;

## Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.61pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CST3

Note: This product is for in vitro research use only